

Recombinant Antibody Expression Services



Fast turnaround



High expression yield



High throughput



High delivery rate



High expression titer

To address a wide range of needs and advance biotechnology, Tsingke introduces the "fast recombinant antibody expression service". Utilizing a high-throughput transient transfection antibody expression platform, we achieve rapid delivery of recombinant antibodies in as fast as one week.


Service Details

Service name	Expression system	Volume	Turnaround time
Recombinant antibody expression service · Anti-dsDNA antibodies · Single chain antibody fragment · VHH-His · VHH-FC · Fab · ...	HEK293/CHO	10 mL	1~2 weeks
		20 mL	
		30 mL	
		40 mL	
		50 mL	
		100 mL	
		200 mL	
		500 mL	
		1000 mL	
		>1 L	

Deliverables

Optional		
· Purified antibodies · COA report(electronic)		· Antibody gene in commercial vectors
Delivery standard	· SDS-PAGE&SEC-HPLC >95%	· Endotoxin level<1 EU/mg

Workflow



Case Study

· Case1




Figure 1. Anti-dsDNA antibody service order

· Case2




Figure 2. VHH-FC antibody service order



Contact us



EN www.tsingke.com | CN www.tsingke.com.cn



market@tsingke.com.cn | info@tsingke.com.cn



Tsingke



tsingkebio



Tsingke Biotech

111240129